

PLASMA STAR

Plasma Xenon Arc Curing Light System

SP-2000

- A New Technological Breakthrough with Style!
- A Powerful Innovation with a Modern Edge for Today's Best Quality Professionals!

- ❖ Digital Functionality Control.
- ❖ Supreme world-class Craftsmanship.
- ❖ High-Energy Efficiency. A time Saver.
- ❖ Great Durability with avant-garde design.



Product Features

- (1) The *PlasmaStar* can produce high-energy plasma arc within spectrum range of **430NM to 470NM** giving the most efficient, yet safely high working power on the patients.
- (2) Innovative **condensing-spotlight** device effectively aggregates light sources creating output beam at over **1800mW/cm²**, yet leaving out unnecessary **infrared and radiant sources to prevent overheat**.
- (3) Original fabrication with authentic Euro-American internal and peripheral components by certified world-class craftsmanship and strict quality control; all under a **3-year warranty**.
- (4) Built-in precision memory device individualizes the suitable operational specifications for different filling and bonding materials. No resetting necessary for each work session.
- (5) Intelligent **LCD backlight design** on the control panel provides much clearer readability, allowing accurate adjustments.
- (6) Internal **spring-switch devising** on the hand-piece gives maximum functionality control without conventional foot pedal.



INSTRUCTION FOR USE

1. SYMBOL

- ⌚ **TIME SET:** It is used to select a cycle of 1, 3, 6, or 10 seconds of curing mode, by depressing the button. The selected cycle will be shown on the LCD display.
- 📶 **LIGHT INTENSITY:** Select light intensity of 50, 80, 100% at normal mode.
- ⚠️ **MODE SETTING:** There are normal mode (NM), ramp curing (RA), 2 Step Curing (SC) and bleaching (BA).
 - NM (normal mode): Fixed power at preset intensity at preset duration.
 - RA (ramp curing): Start at half power then gradually reaches full power for a 6 second duration at the end.
 - SC (2 step curing): 3 seconds half power followed by 3 seconds full power.
 - BL (Bleaching): 3 seconds full power.
- ⊖ **OFF**
- ⊕ **ON**
- ~ **ALTERNATING CURRENT.**

2. SPECIFICATIONS AND EQUIPMENT CLASS

This curing light has been manufactured with respect to electric shock, fire and mechanical hazards only in accordance with UL 2601-1 / can / esa C22.3 no.601-1.

AC supplies connection: 100 to 240 Volts AC (50/60 Hz)

a. Power input: Max. 180 watt

b. Equipment class: Class II (IEC601-1)

c. Protection Against Electric Shock: Type BF (IEC601-1)

d. Protection Against Entry of Liquids: None

e. Use Environment: This equipment is not suitable to use in the presence of a flammable anesthetic mixture with air or nitrous oxide.

f. Operation: At normal mode, there are 4 timed curing cycles of 1, 3, 6, and 10 seconds curing and 50, 80, and 100% intensity output mode.

Fuses (2 Per unit): T5AL250V (110VAC input)

T2.5AL250V (230VAC input)

g. Output wavelength range: 410-500nm (nanometers)

h. Output light intensity: Approx. 1800mw/cm²

i. Standard fiber-optic light guide: 8 mm light liquid guide for general use

j. Overall unit dimensions: Height 10 cm;

Length 35.5 cm;

Width 26.5 cm

k. Unit weight: light-guide 220 g

Control box 4.9 kg